



News – For Immediate Release

Date: February 15, 2016

COX Sterilizer Features Complete 12-Minute Cycles

The COX Rapid Heat Model 6000 provides FDA-Cleared instrument sterilization of wrapped instruments in 12 minutes (6 minutes unwrapped) with no drying cycle.

This HVHA (High-Velocity High Temperature) sterilization system does not require the use of any water, steam or toxic chemicals. Drying cycles are eliminated in HVHA systems and with the absence of water, steam and chemicals there is no corrosion, pitting or dulling of instruments.

The COX Rapid Heat Sterilizer has been recently proven by an independent 3rd party to average 84% less energy per sterilization cycle than competitive steam sterilizers, and can sterilize in double-decker cassettes over 200 instruments per hour. The fast turnaround of the COX Rapid Heat Sterilizer eliminates the temptation of medical staff to bypass or shorten cycles that can jeopardize instrument sterility when instruments are still wet.

In addition to rapid dry heat sterilizers, the company offers a full complement of instrument sterilization pouches, instrument organizers, cassettes, holding racks, spore test kits, incubators and indicator strips.

For more information contact: Diane Baker, dbaker@cpac.com - Toll Free 800-828-6011, extension 1372.
