



News – For December 2022 Release

Maintenance-Free ProVET Sterilizer Doesn't Use Water

RapidHeat™ Sterilization Technology is capturing the attention of the Veterinary Industry with the introduction of a new line of “RH-ProVET” Sterilizers with cycles specifically designed to process veterinary instruments and surgical packs.

Versatile and easy to operate, RapidHeat ProVET Sterilizers are not in the category of autoclave pressure vessels that require water to generate steam. They are an advanced “High-Velocity Hot Air” technology that is 85% more energy efficient than conventional steam sterilizers.

Without water and steam, there are no emissions or residue in the operation of the ProVET sterilizer. It doesn't require an exhaust and drying cycle to maintain load sterility and prevent instrument corrosion. There are no wet loads at the end of a ProVET cycle.

Highly efficient, the RapidHeat ProVET Sterilizer can process instruments “Start-to-Finish” in just 15 minutes. With its spacious rectangular design, the ProVET-11 tabletop sterilizer with its 4 large trays can process as many as 4 large or 8 small surgical packs or 160 pouched instruments in a single cycle.

With its simplicity of touchscreen computerized operation and with no complicated steam and pressure generating components the ProVET sterilizer requires very little or no maintenance.

For more information: visit www.cpac.com, call 800-828-6011, ext. 1372 or E-mail dbaker@cpac.com.